

VX90

User Manual

8290-0000-ENG/V12_EN/2023-01-18



Contents

I - Introduction.....	5
II - Safety.....	7
II.1 - Electricity.....	8
II.2 - Transport, Storage, and Handling.....	8
II.3 - Precautions During Use.....	9
II.4 - Symbols.....	9
III - Equipment and Installation.....	10
III.1 - List of Equipment Supplied.....	11
III.2 - Description of the Device.....	11
2.a) User Side.....	11
2.b) Patient Side.....	12
III.3 - Installation Procedures.....	13
3.a) Unpacking the Unit.....	13
3.b) Electrical Connection.....	14
3.c) Loading Paper into the Printer.....	14
III.4 - Turning the Unit On and Off.....	14
IV - Overview of the Software.....	15
IV.1 - Home Screen.....	16
IV.2 - Keratometry peripheral screen.....	17
IV.3 - Retro-illumination screen.....	19
IV.4 - Corneal Measure.....	20
IV.5 - Data Menu.....	21
V - Performing a Measurement.....	26
V.1 - Preparing to Run a Measurement.....	27
V.2 - Managing the Measurement Procedure.....	27
V.3 - Managing Keratometry Peripheral Measurement.....	27
V.4 - Managing the Retro-Illumination Image.....	29
VI - Configuring the Unit.....	30
VI.1 - Accessing the Configuration Settings.....	31
VI.2 - Overview of the Configuration Screen.....	31
2.a) Main Screen and Refraction Settings.....	31
2.b) Keratometry Settings.....	32
2.c) Device Settings.....	33
2.d) Print Settings.....	34
2.e) Export Settings.....	35
VII - What Should I Do If...?.....	37
VIII - Maintenance.....	39
VIII.1 - Cleaning the Unit.....	40
VIII.2 - Replacing the Printer Paper.....	40
VIII.3 - Replacing Fusible.....	40

VIII.4 - Model eye measurement	41
IX - Appendices	42
IX.1 - Technical Specifications	43
IX.2 - Conformity to Directives and Standards.....	44
2.a) Electromagnetic Emissions	44
2.b) Electromagnetic Immunity	45
2.c) Waste Electrical and Electronic Equipment (WEEE) Directive.....	46
2.d) Manufacturer	47
IX.3 - Contacts Information.....	48

I - Introduction

Product Description:

Thank you for purchasing the Visionix Auto Refractometer VX 90.

The VX 90 is an autorefractometer and a keratometer,

The VX 90 is measuring the objective refraction of the eye giving sphere, cylinder and axis, using an Infra-red ring projection on the Retina.

The keratometer measurement of radius of cornea is done in the same time with another Infra red ring, by measuring the reflection of this ring on the cornea. 4 external fixation dot allow the user to perform peripheral keratometry.

The Machine is showing the intra ocular lens potential opacity with the retro illumination mode.

The practician could also measure the pupil diameter by moving 2 cursors on the screen, the pupil diameter is then displayed on the screen.

A moving fixation target with auto fogging allows the patient to fully desaccomodate in order to get the patient refraction in the best conditions.

Intended Use : The VX 90 is an ophthalmic diagnostic device, autorefractor, keratometer, intended for

- Measuring the refraction of the eye, giving sphere, cylinder and axis
- Measuring the radius of the cornea, central and peripheral keratometry.
- Measuring the pupillary diameter
- Retro illumination for detection of the internal opacities.

Intended Users : the intended users are Opticians, Optometrists, and Ophthalmologists.

The data given by the VX 90 will be used by a qualified person as starting point to start the subjective refraction to make the final prescription. The VX 90 only is not suitable to give the data for the final prescription.

Manufacturer information:



Luneau Technology Operations
2 rue Roger Bonnet
27340 Pont de l'Arche
France

Warning

This document contains confidential information that is the property of Luneau Technology Operations Any use, reproduction or divulging of this material, in part or in whole, is strictly forbidden. This document is provided for the exclusive use of Luneau Technology Operations employees and other authorized users.

The content of this user guide may be modified without warning. The images are not contractual. Every reasonable effort has been made to ensure that its content is accurate. For further information please contact a Luneau Technology Operations representative.

Copyright ©2017 Luneau Technology Operations All rights reserved.

II - Safety

Luneau Technology Operations provides sufficient information to ensure patient safety, avoid system malfunctions, and prevent incorrect readings.

Luneau Technology Operations declines all responsibility for injury to patients or damage to equipment due to ignorance of its safety instructions or in the event that they are not followed.

The safety information appears in the form of warnings and alert messages.

Important!

- Never attempt to disassemble or reassemble the equipment. There are no user serviceable parts in the device.
- Do not modify the equipment in any way.
- Repairs and maintenance must be carried out only by qualified service personnel.
- Operators and patients should keep hands and body clear from moving parts on the device.

II.1 - Electricity

Important!

- To avoid risk of electric shock or bodily injury, do not handle the electrical plugs with wet hands.
- To avoid risk of electric shock or fire, make sure the Vx90's power cord is not damaged before plugging it into an electrical outlet.
- To avoid risk of electric shock, the power cord should be fully inserted in an outlet equipped with a protective ground connection.
- The socket must be easily accessible after the installation of the device
- When connecting external devices such as a screen or printer, ensure that the device conforms to IEC 60950-1 Information Technology Equipment - Safety.

II.2 - Transport, Storage, and Handling

Important!

- Transport the VX90 in its specially designed case.
- Block the VX90 unit with the specially designed screw, below the unit










- Make sure the packing is firm and secure.
- Do not subject the VX90 to strong vibrations. Shocks or violent movements can cause malfunctions.

II.3 - Precautions During Use

Important!

- Do not place or use the VX90 in direct sunlight.
- Do not expose the VX90 to excessive dust or humidity.
- Do not place the VX90 in a hot air current (e.g. above a heater).
- Do not obstruct the ventilation vents.
- Never place the VX90 close to the following types of equipment which can perturb the reception of commands from the remote control:
 - Halogen lamp (direct or indirect)
 - Fluocompact lamp
 - Equipment emitting infrared radiation (auto-phoropter, automatic tonometer, etc.)
- Keep the screen surface clean. Protect it from dust, fingerprints, and shocks.
- When you switch off the VX90, wait at least 5 seconds before switching it on again.

II.4 - Symbols

Symbol	Description
	Important: consult the documents supplied with the equipment
	Type B applied parts
	The equipment must be returned to the manufacturer for scrapping (see <u>Waste Electrical and Electronic Equipment (WEEE) Directive</u>)
	Manufacturer
	Year of Manufacturing
	Serial Number
	Compliance with Medical Device Directive 93/42/EC amended by Directive 2007/47/CE 0051 : identification number of Notified Body IMQ (I)
Rx only	Caution: Federal (U.S.) Law restricts this device to sale by or on the order of a physician.

III - Equipment and Installation


III.1 - List of Equipment Supplied

The VX90 package contains the following items:

- VX90 unit
- Socket head wrench: to unlock plug screw
- Power supply
- Dust cover
- Chin-rest paper
- Printer paper roll
- User Guide
- Model eye
- Fuse


III.2 - Description of the Device

2.a) User Side



N°	Description
1	Optical-measurement head Contains the hardware used to perform all of the measurements during patient diagnostic exams.
2	On/Off switch Shuts down the machine
3	LCD touch screen
4	Up and Down Chin Rest button
5	Joystick to move up /down/left / right the machine
6	Lock Switch To lock the head
7	Reset and Print Key Reset print function
8	Printer

2.b) Patient Side



N°	Description
1	Headrest The patient should lean their forehead on the headrest during all diagnostic procedures.
2	Chin rest The patient should lean their chin on the chin rest during all diagnostic procedures.

III.3 - Installation Procedures

Site Requirements

The VX90 unit should be placed on a clear table or desktop close to a power outlet. The unit should not expose to direct light on the patient side. Better results will be obtained if the unit is located in a room with limited illumination.

3.a) Unpacking the Unit

To unpack the unit:

- Remove the tape around the box.
- take out the carton with care.



- Remove the Polystyrene box to uncover the VX90 unit, which is packed in a protective plastic bag.
- Take the VX90 out of the box and put it on the table.



- Lift the protective plastic bag to uncover the device.

To Unlock the unit:




To unlock the head, turn this screw below the unit

3.b) Electrical Connection

To connect the unit to an electric outlet:

- Check that the power supply voltage corresponds to that required by the equipment (see the identification label on the back of the unit).
- Insert the power connector under the unit
- Connect the power plug to a wall outlet.



under the Vx90

3.c) Loading Paper into the Printer

To insert a paper roll into the printer:

- Lift the handle in the middle of the paper compartment and pull the cover down.
- If an empty paper roll is in the printer, remove it.
- Insert the new roll with the end of the paper on the top of the roll.
- Feed the paper into the slot at the top of the compartment.
- Push the paper-compartment door closed.

III.4 - Turning the Unit On and Off

To turn the unit on:


- Press the On/Off switch (#2 in the illustration).

To turn the unit off:




- On the unit, press the **On/Off** switch.

IV - Overview of the Software

IV.1 - Home Screen




The **Home** screen is the main menu of the VX90. It gives you access to all of the functional screens of the interface.

N°	Description
1	Alignment To select an automatic vertical alignment or a manual vertical alignment: 
2	Measurement Touch to select an automatic or manual measurement 
3	VD Value Touch this key to switch the VD Value: 
4	Cylinder sign To select the cylinder sign, minus, plus or auto.
5	Print mode Use to print ticket or send data
6	Clear Button Press this button to clear the data
7	Settings Menu Touch this button to access the settings menu


8	Refraction Measure Touch this button to select only refraction measurement
9	Keratometry Measure , second push of the button will call the peripheral K mode proceed the same as the central K on the 4 different fixation targets (nasal/temporal/up/down) Third push of this button will call back central keratometry measurement
10	Refraction Keratometry measure Touch this button to have the refraction and keratometry value
11	Keratometry peripheric measurement To open the keratometry peripheric measurement screen
12	Retro-illumination mode To open the retro-illumination screen
13	Corneal Radius It is the mode to measure the corneal radius Second push on the button will call the Retro illumination mode, to acquire the picture of the eye with the back reflexion from the retina. Third push of the button will call back the Corneal radius measurement
14	Data Menu To open the menu of the data you have already measured

IV.2 - Keratometry peripheral screen



N°	Description
1	Print mode Use to print ticket or send data
2	Clear Button Press this button to clear the data
3	Settings Menu Touch this button to access the settings menu
4	Refraction Measure Touch this button to select only refraction measurement
5	Keratometry Measure , second push of the button will call the peripheral K mode proceed the same as the central K on the 4 different fixation targets (nasal/temporal/up/down) Third push of this button will call back central keratometry measurement
6	Refraction Keratometry measure Touch this button to have the refraction and keratometry value
7	Keratometry peripheric measurement To open the keratometry peripheric measurement screen
8	Retro-illumination mode To open the retro-illumination screen
9	Corneal Radius It is the mode to measure the corneal radius Second push on the button will call the Retro illumination mode, to acquire the picture of the eye with the back reflexion from the retina. Third push of the button will call back the Corneal radius measurement
10	Data Menu To open the menu of the data you have already measured


IV.3 - Retro-illumination screen



N°	Description
1	Increase illumination button Used to increase the illumination of the LED
2	Decrease illumination button Used to decrease the illumination of the LED
3	Retro-illumination mode Touch it to have the retro illumination mode
4	Refraction observation in retro-illumination mode Touch it to observe the refraction in the retro illumination mode
5	Print mode Use to print ticket or send data
6	Clear Button Press this button to clear the data
7	Settings Menu Touch this button to access the settings menu
8	Refraction Measure Touch this button to select only refraction measurement

9	Keratometry Measure , second push of the button will call the peripheral K mode proceed the same as the central K on the 4 different fixation target (nasal/temporal/up/down) Third push of this button will call back central keratometry measurement
10	Refraction Keratometry measure Touch this button to have the refraction and keratometry value
11	Keratometry peripheric measurement To open the keratometry peripheric measurement screen
12	Retro-illumination mode To open the retro-illumination screen
13	Corneal Radius It is the mode to measure the corneal radius Second push on the button will call the Retro illumination mode, to acquire the picture of the eye with the back reflexion from the retina. Third push of the button will call back the Corneal radius measurement
14	Data Menu To open the menu of the data you have already measure

IV.4 - Corneal Measure



N°	Description
1	Diameter 1 Touch to modify diameter 1

2	Diameter 2 Touch to modify diameter 2
3	Left perpendicular bar arrow To move the le left perpendicular bar
4	Right Perpendicular bar arrow To move the right perpendicular bar

IV.5 - Data Menu

S	C	A		S	C	A
-1.00	-1.00	000	1	-1.00	-1.00	000
-2.00	-1.50	90	2	-1.00	-1.00	000
+2.00	0.00	000	3	-1.00	-1.00	000
			4			
			5			
			6			
			7			
			8			
			9			
			10			
			CI			
			3.9			

N°	Description
1	Right Refraction measure Display of the all measure with the average at the bottom of the column
2	Left Refraction measure Display of the all measure with the average at the bottom of the column
3	Return Touch it to return at the main screen
4	Refraction Measure

	To select refraction table
5	Keratometry measure To select the Keratometry table
6	Print Touch this button to print one ticket
7	Export Touch this button to export at an external device
8	Delete Touch this button to clear the data
9	Peripheral Keratometry measure To select the Peripheral Keratometry table
10	Confidence Index Give the measure reliability, 4 steps possible: <ul style="list-style-type: none"> - CI=0 All the measures are false - CI=3 the reliability is low - CI=6 the reliability is medium - CI=9 the reliability is High
11	Exam ID reset To reset it long press (about 15 seconds) on "ID".



N°	Description
1	Right Refraction measure Display of the all measure with the average at the bottom of the column
2	Left Refraction measure Display of the all measure with the average at the bottom of the column

3	Return Touch it to return at the main screen
4	Refraction Measure To select refraction table
5	Keratometry measure To select the Keratometry table
6	Print Touch this button to print one ticket
7	Export Touch this button to export at an external device
8	Delete Touch this button to clear the data

The screenshot shows the software interface for the Visionix VX90. At the top, it displays 'R' and 'L' with eye icons, 'ID 0033', and 'PD 30'. The main area is a grid for measurements. The top row is for the Right eye (R) and the bottom row is for the Left eye (L). Each row has columns for R1, R2, A, POS, H-EC, V-EC, and A-EC. The POS column is highlighted in blue. On the right side, there is a vertical toolbar with buttons for Return (R), Keratometry (K), K-P, Print, Export, and Delete. Red boxes and numbers 1 through 9 are overlaid on the screen to indicate specific UI elements.

N°	Description
1	Right Refraction measure Display of the all K-P measurement
2	Left Refraction measure Display of the all K-P measurement
3	Return Touch it to return at the main screen
4	Refraction Measure To select refraction table
5	Keratometry measure To select the Keratometry table

6	Print Touch this button to print one ticket
7	Export Touch this button to export at an external device
8	Delete Touch this button to clear the data
9	Peripheral Keratometry measure To select the Peripheral Keratometry table


Optique Marceau
Paris
01.45.05.05.05
8426

Head: include header (can be modify in the settings), patient ID, Name, Gender, Date, Hour

NAME: _____ M/F
DATE: 2017-04-28 17:29
[REF] UD: 15.00
Cyl. Sign: [-]

<R>	SPH	CYL	AX
	-0.12	-0.12	54
	-0.12	-0.12	55
	-0.12	-0.00	180


AUG -0.12 -0.12 55



Refraction: Include all the value measured on each eyes, the average of each eyes, and diagram of each eyes.

<L>	SPH	CYL	AX
	-0.25	-0.00	180
	-0.25	-0.00	180
	-0.25	-0.00	180

AUG -0.25 -0.00 180



[KER] INDEX: 1.335
<R>

R1	R2	AX
7.94	7.94	180
7.94	7.94	180
7.94	7.94	180

	mm	D	AX
R1	7.94	42.37	180
R2	7.94	42.37	90

AUG 7.94 42.37 -0.00 180

Keratometry: Include all measurement on each eyes and average.

<L>

R1	R2	AX
7.98	7.85	89
7.97	7.85	89
7.97	7.85	89

	mm	D	AX
R1	7.97	42.12	89
R2	7.85	42.75	179

AUG 7.91 42.50 -0.62 89
PD= 72mm

Value of the PD and have the ticket footer, which can be modify in the settings

Visionix

UX90

Ticket:

V - Performing a Measurement

You can initiate a new measurement procedure by opening the Measure screen. The **Measure** screen appears at the start up of the machine

Once the **Measure** screen is opened, you can prepare the patient and select the test options, as explained under Preparing to Run a Diagnostic

V.1 - Preparing to Run a Measurement

After you open the Measure screen, there are different operations you should do before you start running the diagnostic:

- Clean the headrest.
- Place a new sheet of protective paper over the chin rest, or clean the chin rest.
- Instruct the patient to sit down, place their chin on the chin rest, and lean their forehead on the headrest.
- If you need to fine-tune the position of the chin rest, select the **Chin Rest** up and down arrows as necessary to move the chin rest up or down.
- If you need to adjust the position of the unit's head, use the Joystick to adjust the position of the head, up/down / left and right
- Instruct the patient to look at the target (balloon). You should see their eye in the image of the lens view in the **Measure** screen.
- Align the patient's pupil with the center of the target (see Centering the Patient's Eye).

Once you have done all of these things, the measurement start automatically if you are set in automatic mode or manually by pushing the joystick button if you are set in manual mode.

V.2 - Managing the Measurement Procedure

The diagnostic procedure is managed in the Measure screen. In this screen, you can select the measurement you want to perform, position the unit head and chin rest so that the patient's eye is centered in the lens, modify default test settings, start the measurement procedure, and monitor its progress.

Centering the Patient's Eye


In the Measure screen, you can see what the lens is aimed at in the lens view. If the pupil is visible, move the joystick until the pupil is centered and focused. If the patient's pupil is not visible at all in the lens view, adjust the position of the unit's head by moving the joystick to the center of the lens view,

V.3 - Managing Keratometry Peripheral Measurement

It is the mode to measure the curvature of part around cornea. Based upon the center of cornea, measure the curvature of part around cornea from the positions of up/down and left/right direction. It is to indicate the relative eccentricity while comparing the curvature of part around cornea with the curvature of corneal center.

1. Keep pushing K-P button on the right side of the screen.
2. Measurement of Corneal Center

The initial measurement position is the corneal center, and it is indicated as CENTER on the Bottom left of screen. The curvature measured in the corneal center is the same with the one measured in KER mode.



3. Measurement of part around cornea

The direction of part around cornea which is measure at present is to be indicated at the bottom left of measurement mode indication. Four boxes are to be indicated in INF, SUP, NAS, TEMP side of Mire ring. Each box indicates the proceeding state of measurement on part around cornea. If there is the measured result


around part of cornea where the box is located, the inside of box is to be full with color: In case of no result, the box is to be indicated as an empty box. The relevant box indicated at the part around cornea which is measured now is to flicker.

Direction of part around cornea:

- Superior (SUP): Upside from corneal center
- Inferior (INF): Downside from corneal center
- Temple (TEM): To the temple of examinee from corneal center
- Nasal (NAS): To the nose of examinee from corneal center


4. Sequence to measure the part around cornea

Measure it following the sequence of TEM -> SUP -> NAS -> INF In case that the measurement in the direction of part around cornea becomes difficult, the direction lamp (guidance LED light) is to radiate in order to draw the examinee's sight around Mire ring actually. After the examiner shall ask the examinee to look at the light of direction lamp, then he or she can perform the measurement by adjusting the focus of Mire ring.



V.4 - Managing the Retro-Illumination Image.

1. It is to examine the crystalline lens of patient who has the severe symptom of cataract or undergoes it examine the degree of opacity of crystalline lens with the shape of light reflected from retina while changing intensity of light shed on the eye.



1. Adjustment of brightness of the LED
 - Adjust the brightness of the LED clicking on button + and button -.
2. Observation of retro illumination image
 - Place the VX90 in front of the eye (right or left).
 - Use the joystick to center and focus properly.
 - In order to protect the patient's eyes, avoid examining the eyes over 30 seconds.
3. Acquisition of the image
 - Click on the joystick to save an image. It will be saved in the selected box of the current eye. As the live mode of the camera continues after each acquisition, it is possible to take several images consecutively.
 - If more than 3 images are acquired for a same eye, the 1st one will be erased and will be replaced by the 2nd one (and the 3rd one becomes the 2nd one, and the new one becomes the 3rd one).
4. Review an image
 - Click on the wanted image.
 - To come back on live camera, click on the joystick or click on an empty box image of retro illumination.
5. Delete an image
 - Click on the wanted image.
 - Click on clear button. If there is an image after the deleted one, the images will be translated to avoid empty box image.

VI - Configuring the Unit

VI.1 - Accessing the Configuration Settings

To open the Settings Menu:

- Select the **Settings Menu** button.

VI.2 - Overview of the Configuration Screen

The **Configuration** screen gives you access to the VX90's settings and system maintenance functions.


2.a) Main Screen and Refraction Settings



N°	Description
1	Step To choose the step of the measurement, two choice possible 0.12 Diopter or 0.25 Diopter
2	Vertex Distance Touch this button to modify the step of the VD, four choices possible (0, 12, 13.75, 15)
3	Cylinder Notation To choose the sign of the cylinder, plus, minus or automatic
4	# Measures To choose the numbers of measure in one diagnostic (for each eye and for each value)
5	Fogging To choose, if the Fogging is doing at the first measure or all measure when you have selected # Measures > 0
6	Return To return at the main screen
7	Refraction settings To set the taking of the refraction measurement
8	Keratometry Settings To set the taking of the Keratometry measurement


9	Device Settings To change the settings of the device
10	Print Settings To set the printing
11	Export Settings To set all information to export data to an external device

2.b) Keratometry Settings





N°	Description
1	Unit To choose if you want the keratometry in Diopter or in Millimeters
2	Step To modify the step of the measure (0.12 Diopter Or 0.25 Diopter)
3	Index To choose the step of measurement (1.332, 1.336, 1.3375)
4	# Measures To choose the numbers of measure in one diagnostic (for each eye and for each value)
5	Return To return at the main screen
6	Refraction settings To set the taking of the refraction measurement
7	Keratometry Settings To set the taking of the Keratometry measurement
8	Device Settings To change the settings of the device
9	Print Settings To set the printing
10	Export Settings To set all information to export data to an external device

2.c) Device Settings



N°	Description
1	Language Touch the language you want to select
2	Screen Sleep You can select the time before Screen sleeping
3	Beep To choose if you want sound signal when you touch the screen
4	Date and time To select the date and the time of the device
5	Return To return at the main screen
6	Refraction settings To set the taking of the refraction measurement
7	Keratometry Settings To set the taking of the Keratometry measurement
8	Device Settings To change the settings of the device
9	Print Settings To set the printing
10	Export Settings To set all information to export data to an external device


2.d) Print Settings




N°	Description
1	Print To enable or disable the printing of one ticket, or set as automatic
2	R Print To enable or disable the refraction part print in the ticket and choose the type
3	K Print To enable or disable the Keratometry part print in the ticket and choose the type
4	Eye Print To enable or disable the diagram of the eye on the ticket
5	Access to the second page

	Touch this button to access to the second page
6	Access to the first page Touch this button to return at the first page
7	PD Print To enable or disable the print of the PD
8	Ticket Header To choose what you want to write on your own head ticket
9	Ticket Footer To choose what you want to write on your own foot ticket
10	Return To return at the main screen
11	Refraction settings To set the taking of the refraction measurement
12	Keratometry Settings To set the taking of the Keratometry measurement
13	Device Settings To change the settings of the device
14	Print Settings To set the printing
15	Export Settings To set all information to export data to an external device

2.e) Export Settings





N°	Description
1	Export To enable or disable the export data
2	RS 232 to select the protocol of the communication
3	LAN To enable the Wi-Fi connection or disable it
4	VX90 ID Touch "+" "-" to modify the Id of the VX90
5	Access to the second page Touch this button to access to the second page
6	Access to the first page Touch this button to access to the first page
7	WiFi To select the name of the Wifi you want to connect
8	Password To select the password of the wifi you want to connect
11	Return To return at the main screen
12	Refraction settings To set the taking of the refraction measurement
13	Keratometry Settings To set the taking of the Keratometry measurement
14	Device Settings To change the settings of the device
15	Print Settings To set the printing
16	Export Settings To set all information to export data to an external device

VII - What Should I Do If...?

You may encounter some of the following common problems while working with the VX90. If you do, try following the suggested solutions listed below. If the problem persists and the suggested solution does not remedy it, contact a qualified service representative or your local distributor.

Problem	Possible Causes	Suggested Solutions
Images are not properly centered	<ul style="list-style-type: none"> - External lights or windows near the unit - Large patches of light near the unit 	<ul style="list-style-type: none"> - Ensure that there no bright lights or windows near the unit or the patient. - Ensure that there are no large irregular patches of light on the patient or on the unit.
Measured pupil size is different from the expected or actual pupil size	<ul style="list-style-type: none"> - External lights or windows near the unit - Large patches of light near the unit 	<ul style="list-style-type: none"> - Ensure that there no bright lights or windows near the unit or the patient. - Ensure that there are no large irregular patches of light on the patient or on the unit.
Black screen	<ul style="list-style-type: none"> - unit is in stand by mode - unit is switched off - Fuses are blown 	<ul style="list-style-type: none"> - Push on any button to exit the stand by mode - Check the main power switch is on pos 1 - Check the fuses on the main power socket
Tickets do not print	<ul style="list-style-type: none"> - Paper jam - Out of paper 	<ul style="list-style-type: none"> - Open the print paper door and release the paper. - Replace the paper roll.

VIII - Maintenance

This chapter explains how to clean the VX90 unit and how to perform routine maintenance tasks.

VIII.1 - Cleaning the Unit

Important!

Before cleaning the unit, turn it off and unplug it from the electric outlet.

To clean the plastic surface of the VX90, dampen a cloth with a commercial, non-abrasive cleaner and gently wipe the top, bottom, and front surfaces.

CAUTION: Do not spray or pour any liquid directly on the device.

CAUTION: Do not use caustic or abrasive cleaners.

Chin Rest

Cleaning of the chin rest with a soft cloth and alcohol is strongly recommended after each patient. Chin rest paper coverings are supplied with the device. Use the paper to cover the chin rest after each patient use.

Headrest

Regular cleaning of the forehead rest with a soft cloth and alcohol is strongly recommended. The forehead rest is the only part of the machine that comes into contact with the patient.

VIII.2 - Replacing the Printer Paper

The printer is mounted on the right side of the VX90 unit. If the paper roll is used up, the **Printer** LED indicator blinks.

To insert a paper roll into the printer:

- Lift the handle in the middle of the paper compartment and pull the cover down.
- If an empty paper roll is in the printer, remove it.
- Insert the new roll with the end of the paper at the top of the roll.
- Feed the paper into the slot at the top of the compartment.
- Push the paper-compartment door closed.

VIII.3 - Replacing Fusible

Specification : F1AL250V fuse

Position:



VIII.4 - Model eye measurement

The model eye should be measured and the accuracy checked at regular intervals. To set up the model eye, insert the guide groove of the model eye to the chinrest tissue pin.

Set the display step of spherical/cylindrical to 0.12D and perform measurement.

2. Position the model eye

Take down the chin rest paper, alignment the holes of the model eye base and chin rest pad and then insert the chin rest tissue fixed pins.

3. Measuring position adjustment and focus

When the model eye is positioned, move the joystick back and forth to focus until there are eight luminous points appear.

Carefully move the joystick up and down, left and right, until there is a red centering mark appearing.

4. Measuring

As mentioned above, adjust the position of model eye and focus, then press the measure button.



If the measurement result is not within the tolerances given on the model eye, call your dealer.

IX - Appendices

IX.1 - Technical Specifications

Device Specifications

Target Fixation	Auto fog system
Measurement data display	7" TFT color LCD with touch screen (800x480)
Measurement data recording	Built in thermal printer
External output terminal	RS232/WIFI
Power source/power consumption	100V-240V AC, 50/60 Hz, 50VA
Fuses	F1AL250V
Size (L×W×H)	500x288x480 mm
Weight	14Kg
Travel	Back and forth: 40 mm sides: 90 mm Up & Down : 30 mm
Chinrest travel	+/- 30 mm

Operating conditions

Temperature	+10°C to +35°C
Hygrometry	30% to 90%
Air pressure	800 hPA- 1060 hPA

Storage conditions

Temperature	-10°C to + 55°C
Hygrometry	10% to 95%
Air pressure	700 hPA- 1060 hPA

Transport conditions

Temperature	-40°C to + 70°C
Hygrometry	10% to 95%
Air pressure	500 hPA- 1060 hPA

Measurement specifications

Power Range	-20 D to +20 D (step 0.12D/0.25D)
Astigmatism range	0 to 10D (step 0.12D/0.25D)
Axis range	0° to 180°
Axis step	1°
Pupil diameter measured	2 mm minimum (step 0.1 mm)
Retro illumination	IR light
PD Measurement	85 mm measuring range max, 1 mm display unit

Keratometry Specifications

Corneal radius range	5 mm to 10 mm step (0.01 mm)
Corneal refractory power	33D to 67D (n=1.3375) (step 0.12D/0.25D)
Corneal astigmatic power	0 to 12 D (step 0.12D)
Axis range	0° to 180°
Axis step	1°
Measurement diameter	2 mm to 12 mm
Peripheral K	4 external fixation target (nasal/temporal/superior/inferior)

IX.2 - Conformity to Directives and Standards

The VX90 complies with the Medical Device Directive 93/42/EC, as amended by Directive 2007/47/CE and is in Class Im (rationale: rule 12 active devices intended for diagnosis). CE0051 and with the RED Directive 2014/53/EU.

Lifetime of the product: 7 years

1st CE marking: 2019-06

2.a) Electromagnetic Emissions

The VX90 is intended for use in the electromagnetic environment specified below. The customer or the user of the VX90 should assure that it is used in such an environment.


Emission test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11	Group 1	The VX90 uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The VX90 is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purpose.
Harmonic emissions CEI 61000-3-2	Class A	
Voltage fluctuations / flicker emissions CEI 61000-3-3	Complies	

2.b) Electromagnetic Immunity

The VX90 is intended for use in the electromagnetic environment specified below. The customer or the user of the VX90 should assure that it is used in such an environment.

Immunity test	Test level IEC 60601	Compliance level	Electromagnetic environment - guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±(2,4,6) kV contact ±(2,4,8) kV air	Complies	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC 61000-4-4	±2 kV for power supply lines ±1 kV for input/output lines	Complies	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	±1 kV differential mode ±2 kV common mode	Complies	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5 % UT (>95 % dip in UT) for 0,5 cycle 40 % UT (60% dip in UT) for 5 cycles 70 % UT (30 % dip in UT) for 25 cycles <5 % UT (>95 % dip in UT) for 0.5 sec	Complies	Mains power quality should that of a typical commercial or hospital environment. If the user of the VX90 requires continued operation during power mains interruptions, it is recommended that the VX90 be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 hertz) magnetic field IEC 61000-4-8	3 A/m	Complies	Power magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

Note: UT is the a.c. mains voltage prior to application of the test level.

Immunity test	Test level IEC 60601	Compliance level	Electromagnetic environment - guidance
Conducted RF CEI 61000-4-6 Radiated RF CEI 61000-4-3	3 Vrms 150 kHz to 80 MHz 10 V/m 80 MHz to 2,5 GHz	3 V 10 V/m	Portable and mobile RF communications equipment should be used no closer to any part of the VX90, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance $d=1,17\sqrt{P}$ $d=1,17\sqrt{P}$ 80 MHz to 800 MHz $d=2,33\sqrt{P}$ 800 MHz to 2,5 GHz where P is the maximum output power rating of the transmitter in watts (W) according to the manufacturer and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey ^a , should be less than the compliance level in each frequency range ^b . Interference may occur in the vicinity of equipment marked with the following symbol: 

Note 1: At 80 MHz and 800 MHz, the higher frequency range applies.

Note 2: These guidelines may not apply in all solutions. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people

^a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the VX90 is used exceeds the applicable RF compliance level above, the VX90 should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the VX90.

^b Over the frequency range 150 kHz à 80 MHz, field strengths should be less than 3 V/m.

2.c) Waste Electrical and Electronic Equipment (WEEE) Directive



This symbol indicates that the equipment incorporates electronic assemblies and other components that are subject to the Waste Electrical and Electronic Equipment Directive which advise that such electronic and electrical devices should not be scrapped as ordinary domestic waste.

To avoid environmental risks or other dangers caused by irresponsible scrapping, this product and all its accessories must be disposed of separately in accordance with the practices indicated in the WEEE Directive for EU member countries and in local regulations for other countries. For further information on disposal of this product, please contact your local dealer or the manufacturer.

2.d) Manufacturer



LUNEAU TECHNOLOGY OPERATIONS

2 rue Roger Bonnet
27340 Pont De L'Arche
France

IX.3 - Contacts Information

International

VISIONIX INTERNATIONAL SAS
2 Rue Roger Bonnet, 27340 Pont-de-l'Arche - France
Tél. + 33 232 989 132 - Fax + 33 235 020 294
contact@visionix.com
www.visionix.com

Deutschland

VISIONIX DEUTSCHLAND GMBH
An der Pönt 62 - 40885 Ratingen - Deutschland
Tel. +49 (0) 2102-482770 - Fax +49 (0) 2102-48277 77
contact-de@visionix.com
www.visionix.com

España

VISIONIX ESPAÑA, S.A.
Calle Corominas 7 Planta 4 - 08902 Hospitalet de Llobregat - España
Tel. 93 298 07 37
contact-es@visionix.com
www.visionix.com

France

LUNEAU SAS
2 rue Roger Bonnet, 27340 Pont-de-l'Arche - France
Tél. 02 32 98 91 32 - Fax 02 35 02 02 94
contact-fr@visionix.com
www.visionix.com

Italia

VISIONIX ITALIA SRL
Via dei Pioppi 18 - 20024 Garbagnate M.se -Mi-
Tel 02.55413251/221 - Fax 02.55413243
contact-it@visionix.com
www.visionix.com

Portugal

VISIONIX PORTUGAL UNIPessoal LDA
Av. Eng. Duarte Pacheco emp. Amoreiras torre2 - 13ª 1099-042 Lisboa - Portugal
Tel. 214 170 225
contact-pt@visionix.com
www.visionix.com

USA

VISIONIX US
160 Eisenhower Lane North, Lombard, IL 60148
Tel: US: +1 (800) 729-1959
Canada: +1 (905) 760-2420
contact.us@visionix.com
www.visionix.com/us



INNOVATION TO UNLOCK YOUR POTENTIAL

LUNEAU TECHNOLOGY SAS

2 Rue Roger Bonnet, 27340 Pont-de-l'Arche - France

Tél. + 33 232 989 132 - Fax + 33 235 020 294

contact@visionix.com

www.visionix.com